



Image analysis tools: software and libraries used in plant trait measurement

Guillaume Lobet

Université of Liège, PhytoSYSTEMS



@guillaumelobet



Plant Systems Biology
A VIB-UGENT DEPARTMENT





1. What is an image ?
2. What is image analysis ?
3. What is ImageJ ?
4. What are macros and plugins ?




1. What is an image ?

2. What is image analysis ?

3. What is ImageJ ?


4. What are macros and plugins ?

What is an image?




Matrix of values

An image is a matrix of values



		152	230	20	65	98	75	45	25	22	36
	152	230	20	65	98	75	45	25	22	36	65
152	230	20	65	98	75	45	25	22	36	65	36
123	152	148	196	157	249	255	0	0	65	36	65
152	230	20	65	98	75	45	25	22	36	65	36
123	152	148	196	157	249	255	0	0	65	36	65
152	230	20	65	98	75	45	25	22	36	65	36
123	152	148	196	157	249	255	0	0	65	36	65
152	230	20	65	98	75	45	25	22	36	65	36
123	152	148	196	157	249	255	0	0	65	36	65
152	230	20	65	98	75	45	25	22	36	65	
123	152	148	196	157	249	255	0	0	65		

RGB



152	230	20	65	98	75	45	25	22	36
123	152	148	196	157	249	255	0	0	65
152	230	20	65	98	75	45	25	22	36
123	152	148	196	157	249	255	0	0	65
152	230	20	65	98	75	45	25	22	36
123	152	148	196	157	249	255	0	0	65
152	230	20	65	98	75	45	25	22	36
123	152	148	196	157	249	255	0	0	65
152	230	20	65	98	75	45	25	22	36
123	152	148	196	157	249	255	0	0	65
152	230	20	65	98	75	45	25	22	36
123	152	148	196	157	249	255	0	0	65

Greyscale

8-bit
integer [0-255]

32-bit
real values



1. What is an image ?

2. What is image analysis ?

3. What is ImageJ ?


4. What are macros and plugins ?

What is image analysis?





Image analysis **IS NOT** image manipulation
ImageJ **IS NOT** Photoshop

Image analysis is the extraction of information from images



21% of the
leaf surface
is infected




710 cells
average size = 15 px²

Main advantages of automated image analysis


Removes human
appreciation

+

Automation of
processes




area (pixels)




Automated
image
analysis

Phenotyping steps




Types of images

TIFF



530K

JPEG



18K



Minervini et al. (2015). The significance of image compression in plant phenotyping applications. Functional Plant Biology, 1–18


Image analysis basics: Image scale

Principle

Link between pixel
and physical size

DPI

Dots Per Inch
Pixels Per 2.54 cm



Pixels

200

cm

10

scale

**20 px/cm
50 DPI**

Image analysis basics: Histogram

Principle

Distribution of pixel values





Image analysis basics: Thresholding

Principle





Isolate the object from
the rest of the image




Object

Background

Image analysis basics: Thresholding




255



0

Image analysis basics: Thresholding



Be careful with thresholding


use a **fixed** threshold value

or

use an **algorithm**

But use **always** the same

Image analysis basics: Thresholding



Different fixed values





Image analysis basics: Thresholding



Different algorithms




Image analysis basics: **Skeletons**




Image analysis basics: Distance Map





1. What is an image ?

2. What is image analysis ?

3. What is ImageJ ?


4. What are macros and plugins ?

ImageJ

- Open source
- Developed at the NIH
- Created in 1986
- Plugin and macro
- Current version: 1.49

<http://rsb.info.nih.gov>

<http://fiji.sc/>



ImageJ menu

File Basic file operations (opening, saving, creating new images).

Edit Editing and drawing operations as well as global settings.

Image Conversion and modification of images including geometric transformations.


Process Image processing, including point operations, filters and arithmetic operations.

Analyze Statistical measurements, profile and histogram plotting and other operations related to image analysis.

Plugins Commands for creating, editing and managing add-ons

Exercise 1: Thresholding



1. Open the image `bunny.tiff`
2. Duplicate the RGB image
3. Change the image type to 8-bit
4. Duplicate the 8-bit image
5. Threshold the image
6. Save the thresholded image



bit.ly/embo-phenotyping


Exercise 1: Thresholding

1. File > Open
2. Image > Duplicate
3. Image > Type > 8-bit
4. Image > Duplicate
5. Image > Adjust > Threshold
6. File > Save as




Exercise 2: Counting objects

1. Open the image [blobs.gif](#)
2. Set the image [scale](#) to 300 DPI
3. [Threshold](#) the image
4. Create a [binary](#) image
5. [Separate](#) the objects
6. [Count](#) the objects




Exercise 2: Counting objects

1. File > Open Samples
2. Analyze > Set scale...
3. Image > Adjust > Threshold
4. Process > Binary > Make binary
5. Process > Binary > Watershed
6. Analyze > Analyze particles



Exercise 3: Working with roots

1. Open the image [lupin.jpg](#)
2. [Threshold](#) the image
3. Create a [binary](#) image
4. Estimate the [length](#) of the root system
5. Estimate the [diameters of the roots](#)



bit.ly/embo-phenotyping

Exercise 3: Working with roots

1. File > Open
2. Process > Binary > Make binary
3. Image > Duplicate
4. Select [lupin.jpg](#)
5. Process > Binary > Skeletonize
6. Analyze > Analyze particles
7. Select [lupin-1.jpg](#)
8. Process > Binary > Distance Map
9. Process > Image Calculator
10. [lupin.jpg](#) **AND** [lupin-1.jpg](#)





1. What is an image ?

2. What is image analysis ?

3. What is ImageJ ?

4. What are **macros and plugins** ?

Macros and plugins

Macros

Set of ImageJ commands
Useful for automation
Java-like

Plugins

New commands
More complex image analysis
Java

Creating macros

ImageJ built-in macro recording tool

Plugins > Macros > Record...

Macro manual

<http://rsb.info.nih.gov/ij/developer/macro/macros.html>

Launch the macro

Plugins > Macros > Run...

My first macro

```
run("Blobs (25K)");  
  
run("Set Scale...", "distance=300 known=1 pixel=1 unit=cm");  
  
setAutoThreshold("Default");  
  
setThreshold(121, 255);  
  
run("Convert to Mask");  
  
run("Make Binary");  
  
run("Watershed");  
  
run("Analyze Particles...", "size=100-Infinity circularity=0.00-1.00  
show=Nothing summarize");
```

Finish lines with ;
Comment lines with //

A bit more complex

```
setBatchMode(true);


dir=getDirectory("Where are your images");
list=getFileList(dir);
num=list.length;

for(k = 0 ;k < num ; k++){
    open(dir+list[k]);

    run("Set Scale...", "distance=300 known=2.54 pixel=1 unit=cm");
    run("Set Measurements...", "area redirect=None decimal=2");
    setAutoThreshold("Default");
    run("Convert to Mask");
    run("Make Binary");
    run("Watershed");
    run("Analyze Particles...", "size=0.1-Infinity circularity=0.00-1.00
    show=Nothing display summarize");
    close();

}
```

Example of plugin: SmartRoot



More on
Wednesday...



More resources

<http://imagej.nih.gov/ij/>

<http://fiji.sc/>



bit.ly/embo-phenotyping